

Cleaning of charmelle and Q2 brand fabrics



cleaning video

Oil-based stains and other food stains

- Betadine
- Ketchup
- Mustard
- Red wine
- Chocolate

Oil-based stains can be removed with conventional detergents (Persil, washing-up liquid etc.).

- Fresh stains should be removed with a clean cloth. / Dried-on stains should be soaked prior to cleaning.
- Wash off the stain with a warm soapy solution using a sponge.
- Wash off remaining soap (tensides) with plenty of clean water.

Solvent-based stains

- Lipstick
- Nail polish
- Ballpoint pen
- Felt-tip pen
- Shoe polish

Solvent-based stains can be removed with PURE acetone (without additives like gelatine), available from a chemist.

- Carefully apply acetone around the stain.
- Absorb with a soft, clean, white, absorbent cloth as soon as stain dissolves.
- Procedure may have to be repeated.
- Acetone will disappear completely.

Crocking caused by dye-stuff used on other textiles

- Denim
- Other textiles

Discolouring caused by crocking of dyestuff used on other textiles can be removed by using a solution of one part bleach to 10 parts of water.

- Dab discoloured area with the solution.
- Dyestuff particles of the discolouring will be dissolved within 15 – 20 min.
- Absorb the dissolved dyestuff particles with a clean, white cloth.
- Thoroughly rinse out remaining bleach solution.